



Rapid Packaging Strap

ADAPT & OVERCOME TRAINING, LLC

WARNING

FIREFIGHTING, RESCUE, TECHNICAL RESCUE, EMS AND/OR WORKING IN A TACTICAL ENVIRONMENT IS INHERENTLY EXTREMELY DANGEROUS AND MAY RESULT IN SERIOUS INJURY OR DEATH TO THE USER/OPERATOR.

SERIOUS INJURY OR DEATH MAY OCCUR FROM THE IMPROPER USE OF THIS EQUIPMENT

THIS EQUIPMENT HAS BEEN DESIGNED FOR USE BY EXPERIENCED PROFESSIONALS ONLY

USER/OPERATOR MUST THOROUGHLY READ AND UNDERSTAND ALL INSTRUCTIONS CONCERNING BOTH THE CAPABILITIES AND LIMITATIONS OF THIS PRODUCT PRIOR TO USE.

ADAPT & OVERCOME TRAINING, LLC SHALL NOT BE LIABLE FOR ANY DAMAGE, INJURY, OR DEATH RESULTING FROM THE USE OF THE RAPID PACKAGING STRAP. THE RPS USER/OPERATOR ASSUMES ALL LIABILITY, RESPONSIBILITIES OF OWNERSHIP AND RISKS OF USE. KNOWLEDGE IN THE USE OF THE RPS IS A PREREQUISITE OF OWNERSHIP AND IS NOT IMPLIED BY PURCHASE.

THE RAPID PACKAGING STRAP (RPS) IS NOT DESIGNED TO BE A CLASS 1, CLASS 2 OR CLASS 3 HARNESS.

The RPS should not be loaded with more than 300LBS in any configuration.

**MINIMUM BREAKING STRENGTHS (MBS)
System: 1,000 LBS**

Store out of direct sunlight. Direct sunlight will cause breakdown of webbing materials.

RPS Inspection, Cleaning & Retirement

The RPS should be inspected for damage and wear after every use. The RPS is designed for use in harsh environment but should be retired immediately if any of the following are found:

- Rips, holes, abrasions, torn threads, burned, singed or melted webbing.
- Stitching is worn, abraded, or pulling apart
- Exposure to chemicals or corrosive materials including but not limited to: gasoline, diesel fuel, isopropyl alcohol, antifreeze, chlorine bleach, battery acid and/or potentially hazardous exposures from the operational environment.

The RPS is recommended to be retired after 6 years of light use.

If soiled, clean with warm water and mild soap. Allow to drip dry out of direct sunlight. Do not place in a washer machine or heated dryer.



Rapid Packaging Strap

ADAPT & OVERCOME TRAINING, LLC

WARRANTY & RETURN POLICIES

These are the sole and complete warranty and return policies for all products and sets forth the exclusive remedies of the customer. This warranty is given in lieu of all other warranties, either expressed or implied. In no event shall Adapt & Overcome Training, LLC be liable for damages in excess of the product value or for any incidental, special or consequential damages to person or property arising out of the use or inability to use of such product; to the full extent the applicable law may disclaim such damages.

RETURN POLICY

Returns will be handled on a case by case basis. Please contact us at info@adaptovercometraining.com

LIMITED WARRANTY

Adapt & Overcome Training, LLC warrants all its products will be free from defects in construction and workmanship for 90 days from the date of purchase. Adapt & Overcome Training, LLC at its sole option may either repair or replace a defective product at no charge to the customer. Adapt & Overcome Training, LLC is not responsible for damage to the RPS from normal use or misuse by the user and reserves the right to distinguish manufacturing defects from damage caused by the use of the product. A copy of customer purchase receipt and/or invoice may be required for warranty.

CUSTOMER OBLIGATIONS

It is the customer's obligation to fully understand all form, fit and function of the product(s). To obtain warranty service or return services, a customer must complete the warranty Review Request Form included with the product or available online and fax or mail it back to Adapt & Overcome Training, LLC. Once the form is reviewed, a RMA will be issued. The customer must return the product, shipping charges prepaid, along with the warranty request to the provided address.

EXCLUSIONS

The warranty and return policies do not cover any product which has suffered normal wear and tear, defects resulting from accidents, product misuse, product abuse, failure to follow product instructions, improper product maintenance, neglect, tampering with the product, alterations, unauthorized product repairs, fire, flood, chemicals, acts of God, and damage while in transit to Adapt & Overcome Training, LLC. Further, this warranty expressly does not cover any general or consequential damages to person or property other than the product, including but not limited to devices used in conjunction with the product.